

SOLARIS

UNIVERSAL AND POWERFUL BATTERY SYSTEM FOR ORTHOPAEDIC APPLICATIONS

ENGLISH 




GMREIS


Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

SOLARIS is a versatile 3 in 1 battery powered orthopedic surgical equipment: drill, pin driver and sagittal saw. Comfort to the surgeon is provided due: anatomic handle, lightweight, quick and easy attachments coupling reducing the surgical time.



Two drill chuck models: higher rotation speed (a) and higher torque (b).



Two pin driver models: 0.7 to 1.6 mm (a) and 1.8 to 4.0 mm (b).




Sagittal saw with efficient blade coupling mechanism.

Equipment Tray



Autoclavable Components



Non-Autoclavable Components

SOLARIS

CODE	COMPONENT
900-PF-06	Power Drill Handpiece GMReis
900-PF-05	Chuck 4.0mm
900-PF-10	Chuck Key 4.0mm
900-PF-14	Chuck 10.0mm
900-PF-15	Chuck Key 10.0mm
900-PF-01	Pin Driver 0.7 - 1.6mm
900-PF-02	Pin Driver 1.8 - 4.0mm
900-PF-13	Sagital Saw Chuck
900-PF-12	Battery Case
900-PF-08	Transfer Shroud
900-PF-11	Battery
900-PF-07	Charging Platform
900-PF-03	Battery Charger
900-PF-1000	Power Drill Case
711-106	Power Drill Lid
900-PF-1001	Battery Case 300x 200x 100 mm (L x W x H)





GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone .: +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 GMReisBrasil

www.gmreis.com.br