

# STITCH

HIGH RESISTANCE NON-ABSORBIBLE SURGICAL SUTURE AND SUTURE TAPES

ENGLISH 




**GMREIS**

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة


# STITCH

The products of Stitch line - high strength tapes and suture wires are made of UHMWPE polymer that provides high tensile strength and knot resistance, developed specifically for orthopaedic procedures.



USP 2-0, USP 2 and USP 5 surgical suture with adequate strength for each treated region. Different color options facilitate organization during surgical procedures.

All surgical suture and suture tapes are supplied with needle.




Ø3.5 x 15.1 mm Fastlock Anchor with Stitch Tape.



2 mm wide suture tapes provide greater resistance and better performance in contact with soft tissues, ideal for performing the ILA - Internal Ligament Augmentation.


Different color options facilitate organization during surgical procedures.



ILA - Internal Ligament Augmentation with Stitch suture tape and Fastlock anchors for early rehabilitation and early return to sports activities<sup>1,2</sup>.

Excellent biocompatibility and compatibility with metallic and polymeric implants.

High resistance against: breaking, stretching, abrasion, loosening and knots loosening.



<sup>1</sup> Coetzee, J. C., Ellington, J. K., Ronan, J. A., & Stone, R. M. (2018). Functional Results of Open Broström Ankle Ligament Repair augmented With a Suture Tape. *Foot & Ankle International*, 39(3), 304–310. <https://doi.org/10.1177/1071100717742363>

<sup>2</sup> White, W.J., McCollum, G.A. & Calder, J.D.F. Return to sport following acute lateral ligament repair of the ankle in professional athletes. *Knee Surg Sports Traumatol Arthrosc* (2016) 24: 1124. <https://doi.org/10.1007/s00167-015-3815-1>

## UHMWPE Surgical Suture with Taper Point Needle

CODE	MEASURES	COLORS	NEEDLE
333-2-036-AB3817CI	USP 2-0 x 45 cm	Blue and white	3/8 - 17 mm - Cylindrical Tip
333-239-AB1226CI	USP 2 x 90 cm	Blue and white	1/2 - 26 mm - Cylindrical Tip
333-239-PB1226CI	USP 2 x 90 cm	Black and white	1/2 - 26 mm - Cylindrical Tip
333-536-AB1248CO	USP 5 x 90 cm	Blue and white	1/2 - 48 mm - Conventional Tip
333-536-PB1248CO	USP 5 x 90 cm	Black and white	1/2 - 48 mm - Conventional Tip

## UHMWPE Suture Tape with Taper Point Needle

CODE	MEASURES	COLORS	NEEDLE
333-20-AB	2.0 mm x 100 cm	Blue and white	1/2 - 26 mm - Cylindrical Tip
333-20-B	2.0 mm x 100 cm	Black and white	1/2 - 26 mm - Cylindrical Tip

GMReis Stitch surgical suture and suture tapes are supplied sterile for single use, packed to guarantee the integrity of the product.



**HEADQUARTER | GMREIS**

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park  
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone .: +55 (19) 3765 9900  
marketingdigital@gmreis.com.br

   GMReisBrasil

[www.gmreis.com.br](http://www.gmreis.com.br)