


LATARJET

LATARJET FIXATION SYSTEM FOR
TREATMENT OF SHOULDER INSTABILITY

ENGLISH 



GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

2022


LATARJET


Latarjet technique consists of removing a coracoid graft from scapula to fix on anterior border of glenoid, thus enlarging the articular contact surface with the humeral head to reduce shoulder instability and dislocation.

Latarjet GMReis System was developed to provide the surgeon with a safe and effective technique for both graft removal and positioning the fixation in glenoid, using the Hashi technique.

Graft fixation is performed with two cannulated screws developed specifically for Latarjet system, with diameter and lengths suitable for implantation region, available full or partial thread options.

To increase the screw contact surface with the graft, improving compression performance, Latarjet GMReis System offer washers or plate with two specific holes.





LATARJET CANNULATED COMPRESSION SCREW

CODE	SCREW THREAD	LENGTH
324-385-30	total	30.0 mm
324-385-32	total	32.0 mm
324-385-34	total	34.0 mm
324-385-36	total	36.0 mm
324-385-38	total	38.0 mm
324-385-40	total	40.0 mm
324-385-42	total	42.0 mm
324-385-15-30	15 mm	30.0 mm
324-385-15-32	16 mm	32.0 mm
324-385-15-34	17 mm	34.0 mm
324-385-15-36	18 mm	36.0 mm
324-385-15-38	19 mm	38.0 mm
324-385-15-40	20 mm	40.0 mm
324-385-15-42	21 mm	42.0 mm

LATARJET INSTRUMENTAL KIT

CODE	DESCRIPTION
324-15	Latarjet Guide Clamp
324-10	Latarjet Guide
324-12	Screw Positioning Wrench
324-11	Depth Gauge
324-13	Drill Ø2.7 mm x 195.0 mm Cannulated 1.2 mm
324-14	Drill Ø2.7 mm x 265.0 mm Cannulated 1.2 mm
324-16	Drill Ø4.0 x 120.0 mm
324-17	Guide Wire Ø1.0 x 280.0 mm
324-18	Guide Wire Ø1.0 x 210.0 mm
324-ES	Latarjet Case



FLAT SUTURE WASHER

CODE

324-385-30



LATARJET SUPPORT PLATE

CODE

324-01



LATARJET DISPOSABLE KIT

CODE	MODEL
324-CX-500-100	SHORT
324-CX-500-200	LONG



HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone : +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 **GMReisBrasil**

www.gmreis.com.br