

# REDUX COLLINEAR CLAMP

MINIMALLY INVASIVE FRACTURE REDUCTION CLAMP

ENGLISH 



**GMREIS**


Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

2021

# REDUX - COLLINEAR CLAMP

GMReis Collinear clamp were developed to assist the surgeon in reducing and maintaining fragment position minimally invasive during fractures surgeries of: joints, long bones (tibia and femur) and pelvis.



*Minimally invasive reduction of condylar fracture with Percutaneous Redux Clamp.*




*Pelvic Redux Arm for minimally invasive pelvis fractures reduction*

*Two fixing options:  
Threaded Linear Guide  
Smooth Linear Guide*



*Redux Arm Hohmann type for minimally invasive reduction fractures of long bones (tibia and femur)*




**COLLINEAR GUIDE KIT**

CODE	DESCRIPTION
901-100S	Disposable Collinear Guide Kit

**REDUX INSTRUMENTAL KIT**

CODE	DESCRIPTION
232-236	Redux - Grip of Collinear Reduction
232-108	Redux Arm Hohmann type
232-110	Pelvic Redux Arm
232-111	Percutaneous Redux Arm
901-107	Cleaning Thread
901-106	Cleaning Brush



**GMREIS**

*Qualidade para Vida*

*Calidad para Vida Quality for Life الجودة للحياة*

**HEADQUARTER | GMREIS**

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park  
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone .: +55 (19) 3765 9900  
marketingdigital@gmreis.com.br

 GMReisBrasil

[www.gmreis.com.br](http://www.gmreis.com.br)