


LOW PROFILE UNI ELBOW

A NEW CONCEPT FOR INTERNAL ELBOW STABILIZATION

ENGLISH 



GMREIS

Qualidade para Vida


Calidad para Vida Quality for Life الجودة للحياة

2022

UNI ELBOW

Uni Elbow^{GMReis} was developed for the treatment of unstable elbow injuries through internal fixation temporary.

Uni Elbow^{GMReis} provides excellent stabilization, allowing early active movement, reducing the inconveniences and complications associated with external fixation.









Fig.: Elbow in 90° flexion and extension with the Uni Elbow internal stabilizer.

CANNULATED FIXING PIN


CODE	LENGTH
345-12-30-C	30 mm
345-12-35-C	35 mm
345-12-40-C	40 mm
345-12-45-C	45 mm
345-12-50-C	50 mm
345-12-55-C	55 mm
345-12-60-C	60 mm
345-12-65-C	65 mm
345-12-70-C	70 mm



CODE	DESCRIPTION
345-10	Elbow stabilizer

UNI ELBOW INSTRUMENTAL KIT

CODE	DESCRIPTION
345-18	Axis Centering Guide
345-19	Axle Guide
345-20	Wire Guide
329-207	Drill Ø2.5 x 150.0 mm Graduated
345-14	Cannulated Drill Ø2.7 x 70 mm
345-15	Torxdrive T10 Wrench Fixed with Silicone Handle
345-16	Small Depth Gauge
169-154	Depth Gauge 10-60mm
172-19A	Wire Cutter Ø3.2 mm
223-323	Guide Wire Ø1.0 x 150.0 mm Without Thread





HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone : +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 GMReisBrasil

www.gmreis.com.br