

CARPOLUX

OPTICAL CARPAL TUNNEL RELEASE

ENGLISH 



GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

2021


CARPOLUX

The Carpolux^{GMReis} is a surgical device developed for carpal tunnel release surgical procedure, supplied sterile for single use.

Carpolux's built-in lighting system provides better visualization of the treated regional lowering safe and effective technique, reducing the surgical time and soft tissue morbidity.



Fig.: Insulated cutting blade by the atraumatic tip of the Carpolux.



Carpolux lighting activation device



FEATURES AND ADVANTAGES

- Safe and effective technique, with an incision of only 2 cm, it is possible to visualize the distal third of the transverse ligament;
- Ligament lighting;
- Cutting blade isolated from surrounding soft tissues, avoiding unintentional injuries;
- Atraumatic tip to prevent injury to nerves and,
- Procedure can be performed under local anesthesia, without the need for additional equipment..

CODE	DESCRIPTION
265-01	Carpolux - Carpal Tunnel Release System.



HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone .: +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 GMReisBrasil

www.gmreis.com.br